


Clinisys Atlas application installation instructions for Shared Print Service.

- **Pop-Up Blockers should be disabled for the URL used to connect to the Clinisys Atlas application.**
- **This version supports Windows operating systems.**


Clinisys Atlas URL: <https://cellnetixsso.elaborders.com>

To manually download the SPS installer file:


1. Enter the .MSI download URL in your browser window:
<https://testing-cellnetix.elaborders.com/Implementation/SPSApplication.msi>
2. This will prompt you to save the file to your local computer (select a file or save in a default folder, such as “Downloads”).




3. Navigate to the folder where you saved the file and double-click it to begin the installation process. Note, you may need to click on alerts based on your Windows security settings.
4. The Atlas Print Service Setup Wizard will appear:




5. Click “Next” to proceed to the next screen:




6. Click "Install" to begin. This will show the installation screen with a green progress bar:




7. Once the installation is done a new screen will appear. Select the Default printer for SPS to use for printing documents out of the Sunquest Atlas application. Also, you may choose a directory to save PDF files printed via SPS, or you may keep the default directory provided by the application.



8. Click "Next" after making your selections. You will then be presented with the option to perform a test print to a selected printer:




9. Click “Next” to bypass this step and test printing directly from the Sunquest Atlas application. Or, select a printer and click “Test Print” to perform a test print.



10. You will now see the SPS Application icon on your system tray to show that it's running:



11. Right-click on the SPS icon to open the PDF directory, change the default PDF directory or change the default printer:



12. With your default printer selected you should now be able to print documents from the Sunquest Atlas application as needed, when you access a site that is enabled for SPS printing.